

APPLICATION POINTS

Categories of illnesses depending on body part and affected area	Programme name	Applicator type	Place of application
Internal organs	goat, H1, H2	CAVO	affected area
	inflammation of the urinary tract, H2	MANIC, CAVO	lower abdomen
	inflammation of the pancreas, H2	LETTINO, VIERA, PIANO	area between thoracic (chest) and lumbar (lower back) spine
	inflammation of the ovaries and fallopian tubes, H2	MANIC, CAVO	lower abdomen
	stomach (gastric) ulcers, H2	LETTINO, VIERA, PIANO	area around stomach
	inflammation of the middle ear (otitis media), H2	CAVO	head

APPLICATION POINTS

Categories of illnesses depending on body part and affected area	Programme name	Applicator type	Place of application
Musculoskeletal system	joint arthrosis, H1	CAVO	affected joint
	frozen shoulders, H2	CAVO	shoulders
	back and spinal pain, H2	SIECCO, LETTINO	spine
	carpal tunnel syndrome, H2	CAVO, MANIC	wrists, cervical spine
	cervical spine blockage, H2, H4	CAVO	cervical spine
	osteoporosis, H1	CAVO, SIECCO PIANO, LETTINO	affected area
	venous ulcers, H1, H2	CAVO, PIANO LETTINO, VIERA	affected area
	spinal bone spurs, H2, P1, P2	SIECCO, LETTINO	spine
	rheumatoid arthritis, H2	CAVO, SIECCO, PIANO, LETTINO	spine and affected area
	herniated disc, H2, P1, P2	SIECCO, LETTINO	spine
	fractures - complicated healing, H1, H2	CAVO, LETTINO, VIERA	affected area
	bone spurs, H2	CAVO, MANIC	heel and foot
	diabetic foot, H1, H2, H4	CASSO, LETTINO, PIANO	affected area
	limb injuries, H1, H4	CAVO, LETTINO	affected area
	Sciatica, H2	CASSO, LETTINO	spine
Bechterew's disease, H2	CASSO, LETTINO	spine and affected area	

APPLICATION POINTS

Categories of illnesses depending on body part and affected area	Programme name	Applicator type	Place of application
Musculoskeletal system	aseptic bone necrosis, H4	CAVO, LETTINO	affected joint
	Endoprosthesis, H2	CAVO, VIERA	joint
	Lumbago, H2	CASSO, LETTINO, PIANO	spine
Cardiovascular system	venous diseases, H1	CASSO	whole body
	arterial diseases, H1	CASSO	whole body
	high blood pressure (hypertension), H1	LETTINO, PIANO	chest area
	lymph vessel diseases, H2	CASSO	whole body
	promotion of healing after surgery, H2, H4	CAVO, LETTINO, VIERA, CASSO	affected area
	swelling, H1, H2	CAVO, LETTINO, VIERA, CASSO	affected area
	Neurology	Alzheimer's disease, H1, H2	CAVO, LETTINO, CASSO
Parkinson's disease, H1, H3		CAVO, LETTINO, VIERA, CASSO	head, spine
peripheral neuropathy, H4		CAVO, LETTINO, CASSO	head, spine
tinnitus, H1, H2		CAVO	head
multiple sclerosis, H2		CAVO, LETTINO, CASSO	head, spine
glaucoma, H2		CAVO	head

APPLICATION POINTS

Categories of illnesses depending on body part and affected area	Programme name	Applicator type	Place of application	
Neurology	tick-borne encephalitis, H2, H4	CAVO, LETTINO, CASSO	head, spine	
	lyme disease, H2, H4	CAVO, LETTINO, CASSO	head, spine	
	Migraine, H1, H2	CAVO	head, spine	
	damage caused by stroke, H3	CAVO	head, spine	
	optic nerve atrophy, H3, H4	CAVO	head	
	polio, H3	CAVO, LETTINO, CASSO	head, spine	
	Myopathy, H1	LETTINO, CASSO	spine	
	Internal organs	Diabetes, H1	LETTINO PIANO	area around pancreas and liver
		erectile and sexual dysfunction, H1, H2	MANIC, CAVO	genital area
Haemorrhoids, H2		MANIC	anal area	
hypothyroidism (under-active thyroid), H3, H4		CAVO	larynx	
immune system boost, H1, H2		CAVO	area around thymus	
prostate inflammation, enlargement, H4		MANIC, CAVO	lower abdomen	
inflammation of the kidneys, H2		LETTINO, VIERA	lumbar region	
promotion of healing after surgery, H2, H4		CAVO, LETTINO, VIERA, CASSO	affected area	